

# Mezoforte® Duo

## Better penetration of active agent by electroporation and skin tightening by radiofrequency

Several methods are known for the facilitating the absorption active products into the skin. The skin layers which block the entry of active ingredients are removed by mechanical or enzymatic peeling or microdermabrasion. Using electricity may increase the penetration of the products, but entry via needle has proved more efficient recently. However, a significant proportion of clients don't prefer the needle method.

Now Dr Derm Equipment presents the mezophoresis® treatment which utilizes the new needling mesotherapy method. In the course of the new method of entry of active ingredients, the double lipid layer, which represents the barrier of the skin, is rearranged owing to the electric pulses that are in the millisecond range in duration. Such structural change results in the formation of new channels for the active ingredients, thus the skin readily absorbs the active molecules applied onto its surface.

The bipolar RF handpiece is used for wrinkles of the face and other treatments requiring superficial effect. In this case the radio waves flow between the two poles built into the handpiece, thus primarily reaching the epidermis and the upper dermis, providing rejuvenation of the skin surface.

The Mezoforte® Duo unit continuously controls the performance of radiofrequency. Owing to *Computerized Temperature Control of Radiofrequency (CTC<sup>RF</sup>)*, it protects the epidermis from overheating, and thus from the inconveniences associated with some other devices.

The Mezoforte® Duo offers two handpieces and technologies together or separately. As usual with Dr Derm devices, the control of Mezoforte® Duo is fully digitalized, which provides a high level of treatment safety. The device offers appropriate programs for each handpiece, thus the user needs only to select the area to be treated, and the treatment time. The intensity of treatment may be adjusted according to the given problem and the sensitivity of the skin.

The operators of Dr Derm Equipment can modify the new device setting according to their experiences, the radiofrequency treatment is painless and easy to perform.

The mezophoresis® method provides a more intensive effect due to faster and more efficient absorption of the active ingredients applied onto the skin.

Thus, Mezoforte® Duo the device offers the latest technology with the usual stability of value, in excellent Dr Derm quality, with a 3-year warranty.



Superbrands  
2009

Dr Derm Equipment Kft.  
H-1024 Budapest, Fény u. 2. • Telefon: +36-70 770-2664  
drderm@drderm.net • www.drderm.net

DRDERM  
EQUIPMENT